

## iSeries Bed – Semi or Full



### I Series Bed Overview

Two models: full electric (Bed 1300I) and semi electric (Bed 1200I)

Neither model includes bed rails

Half rail (pair): B-Rail120I

Full rail: B-Rail230I

Unique 3-in-1 capability: converts between full electric, semi electric, or low bed

### Unboxing and Assembly

Ships in two boxes: head and foot decks sandwiched together in one; second box contains head and footboards.

Open top, flip box upside down, slide cardboard off the top.

Head and foot sections zip-tied together with foam padding bring scissors or wire cutters.

Tools required: hex T wrench or hex Allen wrenches

Motor arrives detached; one-time attachment takes ~5 minutes. Once the second motor is installed post shipping, dealer never needs to remove it. Should check allen bolts occasionally to make sure they remain secure.

If they purchase configured semi electric, the upgrade kit will contain the second motor and full electric control pendant.

## iSeries Bed – Semi or Full



- Assembly of head and foot deck is flat (no 90-degree tilting needed)
  - Slightly lift decks, insert male into female receptacle on both sides
  - Push together using bar deck, insert cotter pin to lock
- Headboard and footboard labeled with small stickers (head/foot): check before installing.
- Assemble head and foot section decks together first. Step to either end then lift up the end of the deck. Insert both pegs into deck hooks at the bottom by angling head/footboard slightly outward. Once bottom pegs are in place, bring to 90-degree angle and the deck locks into place on both pegs.
- Integrated locks on headboard and footboard: slide outward to lock permanently in place. Lift bed end by deck to demo locks are working.
  - Unique to I series; competitors use toggle bolts or cotter pins.

### **Key Design Features to Highlight.**

**Be sure to point out thickness of steel, smooth welds, and how robust the deck hooks are.**

Dual actuator arms on head and foot sections

Two arms per section vs. one on older beds

Reduces noise, adds longevity, eliminates uneven sag when fully raised

Steel slat deck vs. spring deck

No spring compression needed to join decks; springs wear and dip over time; steel is more rigid and durable.

Slats have lower chance to rip mattress covers.

Hooks relocated to the deck (not headboard/footboard)

Easier transport: stack flat against truck wall or warehouse shelves

Stronger welds as integrated part of deck

Pegs now on inside of headboard/footboard leg for added strength

Four locking rubber-edged casters (not cheap plastic; no mix of locking/unlocking)

All headboards and footboards include high/low posts

Backup hand crank for full electric clipped to footboard

Quiet operation: high-strength motor with plastic grommet on tension rod to insulate noise.

## iSeries Bed – Semi or Full



### **Motor Installation and Controller**

Two cylinders on main motor face outward (same side as controller and second motor plug)

Motor marked “head” and “foot” for correct orientation

Remove slide-over caps, maneuver motor onto bar (raise head or foot section slightly if needed), slap from bottom to seat, replace caps

Plug foot motor phone cord in (one port only), then plug in controller

Controller functions:

- Up/down for full bed height

- Up/down for head section

- Up/down for leg section

- Raise all / lower all simultaneously (fourth function, unique to I series)

Tension rod height settings: two positions, one for semi electric, one for full electric

- Full electric: shorten rod to lowest setting